

校正報告書

校正暨量測實驗室-台北

Report Date: 2024/6/13

	Report No.: E	CR2421960	第1頁	共2頁
	臺-	比醫學大學附設醫院		
		-20℃冷藏櫃		
LIEBHERR	機型 Model	EFE-1102	序號 Serial No.	6001023-000131
GENP-EC-T014(V 1.6)	收件日期 Received Date	2024/6/6	校正日期 Calibration Date	2024/6/6
耿志昌	溫度 C Temperature	(20 ~ 19) °C	相對濕度 % Relative Humidity	(67 - 65) %
110臺北市信義區吳興街252號 11031台北市信義區吳興街252號蔡萬才大樓		校正地點 Location	遊枚	
	LIEBHERR GENP-EC-T014(V 1.6) 耿志昌	基: LIEBHERR 機型 Model GENP-EC-T014(V 1.6) 收件日期 Received Date 耿志昌 溫度 C Temperature 110臺北市信義區吳興街2	臺北醫學大學附設醫院 -20°C 冷藏櫃 LIEBHERR 機型 EFE-1102 Model GENP-EC-T014(V 1.6) 收件日期 2024/6/6 Received Date 耿志昌 溫度 C (20~19) C Temperature 110臺北市信義區吳興街252號	臺北醫學大學附設醫院 -20℃冷藏櫃 LIEBHERR 機型 EFE-1102 序號

儀器名稱 Equipmont	製造廠商	機型 Madal	標準器校正日期
Equipment	Manufacturer	Model	Calibration Date
Thermometer	СНУ	506R	2023/8/2
序號	追溯單位	報告號碼	標準器有效日期
Serial Number	Traceability	Report No.	Due Date
50000511	SGS(TAF 0143)	ECR2329663	2024/8/1

- ◆ 台灣檢驗科技股份有限公司特此聲明本報告書內記載之標準器,依ILAC PEO之規定,可追溯至ILAC MRA 國際實驗室認證聯盟相互采認協議或員 之仍可實驗室,或各國家計量標準機構(NMI),或國際度量衡委員會相互認可協定之機構(CIPX MR4)。或驗證本有物質(CRXS)。 有TAF犯避操结之報告表示符合全國認證基金會(TAF)認證範围之校正項日證書;無TAF認證標誌之報告亦符合本實驗室標準校正作業程序及1S0 [EC 17025之規定。
- SGS Tarwar Ltd hereby declare that traceability follows ILAC P10(ILAC policy on the traceability of measurement results) All standards are directly traceable to TAF recognized lab (members of the ILAC MRA) or to National Metrology Institutes (NMI) or to other international standards (members of the CIPM MRA) or cartified reference materials (CRMs) The report with "TAF" accredited symbol represents the certificate of calibration items conforming to TAF accreditation scope , Writhout "TAF" accredited symbol , the report also complies with the lab s standard calibration operating procedures and ISO/IEC 17025 requirements
- 本校正報告僅對上遞儀器之校正項目有效且未經調整,報告內容不包含校正週期之建議以及符合性判定、除特別僱明外。
- The calibration report is only valid for the instrument mentioned above and has not been adjusted Unless otherwise specified the content of the report does not include the suggestion for calibration interval and the determination of compliance
- 本校正報告部份權製及影本無效。
- To reproduce or copy calibration report in partial is not allowed
- 校正程序名稱: 烘箱 低温箱 實測 校正作業程序書 Procedure used
 - Calibration SOP for Oven (Actual test)

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com.tw/Terms-and-Conditions and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com.tw/Terms-and-Conditions. Attention is drawn to the limitation of liability, and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not experient be a transaction for exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company Any unauthorized alteration, forgery or refer only to the sample(s) tested.

TWE 0918151



校正報告書 (Calibration Report)

RptNo.: ECR2421960

第2頁 共2頁

温度部份:(比對校正)

器示值(℃)

標準值(℃)

器差值 (℃)

-21.1

-21.2

0.1

校正說明:

1,器差值=器示值-標準值

- 2.器示值係指送校正件所顯示或設定之值
- 3.標準值係指工作標準件之輸出值或顯示值
- 4.校正能力係以約95%信賴水準,k=2之相對擴充不確定度表示
- 5.擴充不確定度:0.5℃
- 6.上述擴充不確定度已包含校正件之不確定度評估結果 -- THE END --

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com.tw/Terms-and-Conditions and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com.tw/Terms-and-Conditions. Attention is drawn to the limitation of liability, and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or refer only to the sample(s) tested.