

校正報告書

校正暨量測實驗室-台北

		Report No.: E	CR2421962	第1頁	共2頁
申請者 Applicant		臺:	比醫學大學附設醫院		
機器名稱 Equipment			4℃冰箱	(>h \$\overline{1}{1}\overline{1}{2} B)	
製造廠商 Manufacturer	Thermo Fisher Scientific	機型 Model	FRCR4504A	序號 Serial No.	015597020116121
校正程序 Procedure used	GENP-EC-T014(V 1.6)	收件日期 Received Date	2024/6/6	校正日期 Calibration Date	2024/6/6
校 验者 Operator	耿志昌	温度 [®] C Temperature	(20 ~ 19) °C	相對濕度 % Relative Humidity	(67 ~ 65) %
顧客地址	110畫:	比市信義區吳興街2	52號	校正地點	·
遊校地址	11031台北市	信義區吳興街252號	蔡萬才大樓	Location	
	實驗	(室使用標準器	/ SGS Standar	ds	
		製造廠商			校正日期
Equi	pment Mar	mfacturer	Model	Calibra	tion Date
There	nometer	СНУ	506R	202	3/8/2

序號	追溯單位	報告號碼	標準器有效日期
Serial Number	Traceability	Report No.	Due Date
50000511	SGS(TAF 0143)	ECR2329663	2024/8/1

- ◆ 台灣檢驗科技股份有限公司特此聲明本報告書內記載之標準器,依ILAC P10之規定、可追溯至ILAC MR4 國際實驗室認證聯盟相互承認認填成員 之認可實驗室、或各國家計畫標準機構(YEI),或國際產量衝安員會相互認可論定之機構((IPX MRA),或驗證本考物質((RMS)。 有TAF認證標誌之報告表示符合公園認證基金會(TAF)認證範围之校正項目證書;無TAF認證標誌之報告亦符合本實驗宣標準校正作業程序及ISO/IEC 17025之規定。
- SGS Taiwan Ltd hereby declare that traceability follows ILAC PublicIAC policy on the traceability of measurement results) All standards are directly traceable to TAF recognized lab (members of the ILAC MRA) or to National Metrology Institutes (NMI) or to other international standards (members of the CIPM MRA) or certified reference materials (CRMs). The report with "TAF" accredited symbol represents the certificate of calibration items conforming to TAF accreditation scope , Without "TAF" accredited symbol , the report also complies with the lab si standard calibration operating procedures and ISO/IEC 17025 requirements
- 本校正報告僅對上遠儀器之校正項目有效且未經調整,報告內容不包含校正週期之建議以及符合性判定,除特別聲明外。
- The calibration report is only valid for the instrument mentioned above and has not been adjusted Unless otherwise specified, the content of the report does not include the suggestion for calibration interval and the determination of compliance
- 本校正報告部份複製及影本無效。
- To reproduce or copy calibration report in partial is not allowed
- Procedure used
- 烘箱/低温箱(實測)校正作業程序書 Calibration SQP for Oven (Actual test)

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com.tw/Terms-and-Conditions and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com.tw/Terms-and-Conditions. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Chent's instructions, if any. The Company's sole responsibility is to its Client and this document does not expert partials to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company Any unauthorized alteration, forgery or refer only to the sample(s) tested.



校正報告書 (Calibration Report)

RptNo.: ECR2421962

第2頁 共2頁

温度部份:(比對校正)

器示值(℃)

標準值(°C)

器差值(℃)

4.8

4.8

0.0

校正說明:

- 1.器差值=器示值-標準值
- 2.器示值係指送校正件所顯示或設定之值
- 3.標準值係指工作標準件之輸出值或顯示值
- 4.校正能力係以約95%信賴水準,k=2之相對擴充不確定度表示
- 5. 擴充不確定度: 0.5℃
- 6.上述擴充不確定度已包含校正件之不確定度評估結果 -- THE END --

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com.tw/Terms-and-Conditions and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com.tw/Terms-and-Conditions. Attention is drawn to the limitation of liability, and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not expone a proper a parties to a transaction from exercising all their relations on the content or appearance of this document cannot be reproduced except in full, without prior written approval of the Company Any unauthorized alteration, forgery or refer only to the sample(s) tested. TWE 0918154